CONDUCELL 4USF ARC SENSORS

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Important note

Copyright © 2010 HAMILTON Bonaduz AG, Bonaduz Switzerland. All rights reserved. The reproduction of any part of this document in any form is forbidden without the express written agreement of HAMILTON Bonaduz AG.

Contents of this manual can be modified without previous announcement. Technical modifications reserved. Greatest possible care was used on the correctness of the information in this manual. If errors should be discovered nevertheless, HAMILTON Bonaduz AG is pleased to be informed about it. Regardless of this, HAMILTON Bonaduz AG cannot assume liability for any errors in this manual or for their consequences.

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1. INTRODUCTION	4
2. INTENDED USE	4
3. SAFETY INSTRUCTIONS	4
4. INITIAL OPERATION	5
5. ELECTRICAL CONNECTION	5
ELECTRICAL CONNECTION OF THE 4–20 mA CURRENT INTERFACES	6
EXAMPLES OF CIRCUIT ARRANGEMENT	7
ELECTRICAL CONNECTION FOR THE DIGITAL RS485 INTERFACE	8
EXAMPLE OF A CIRCUIT ARRANGEMENT	9
6. CONFIGURATION AND MONITORING OF THE SENSOR	10
7. PREPARATION FOR THE MEASUREMENT	10
8. REMOVAL OF THE SENSOR	11
9. STERILIZATION, AUTOCLAVING, CIP PROCEDURES	11
10. TESTING AND MAINTENANCE	12
CALIBRATION	12
SELF-DIAGNOSIS FUNCTIONS	13
11. DISPOSAL	14
12. PARTS AND ACCESSORIES	15
13. TECHNICAL DATA	16
BEDIENUNGSANLEITUNG DEUTSCH	17

Operating Instructions CONDUCELL 4USF ARC

1. Introduction

This manual refers to the CONDUCELL 4 USF ARC range of sterilizable conductivity sensors from HAMILTON Bonaduz AG. These sensors are compatible with all other components of the HAMILTON ARC System, including a complete family of intelligent sensors for pH, dissolved oxygen and conductivity measurements in a process control.

Designation	Order number P/N
CONDUCELL 4USF ARC 120	242159
CONDUCELL 4USF ARC 225	242160
CONDUCELL 4USF ARC 325	242161
CONDUCELL 4USF ARC 425	242162

Hamilton ARC sensors are quality products manufactured according to the very latest research findings. Follow the instructions given here to be sure of optimal safety and durability.

IMPORTANT: These instructions must be read, understood and followed by all staff using CONDUCELL 4USF ARC sensors. HAMILTON can assume no responsibility for damage and operational disruption arising from failure to observe these instructions.

2. Intended use

CONDUCELL 4USF ARC sensors were developed for the measurement of the conductivity in liquids. The engineering design based on 4-electrode principle guarantees an excellent linearity.

Main application for CONDUCELL 4USF ARC sensors is the in-line measurement under demanding conditions in biotechnology and the general chemical industry. During development, special attention was paid to an optimum sanitary design. The sensor meets the hygienic standards of EHEDG (Hygienic Equipment Design Criteria). All materials in contact with the solution are FDA approved.

CONDUCELL 4USF ARC sensors provide both standard analog (4–20 mA) and digital Modbus interfaces. These are built into each sensor and are supported directly from the sensor head. They do not require any additional equipment such as amplifier or transmitter.

3. Safety instructions

CONDUCELL 4USF ARC sensors must be used for their intended applications, and in optimum safety and operational conditions. The specifications (such as temperature or pressure) defined in the section entitled 'Technical data' must not be exceeded under any circumstances. Potential hazards can exist if the sensor is not operated correctly or appropriately.

Assembly and maintenance must be performed only by trained personnel.

Make sure that the PG 13.5 thread and the O-ring are not damaged when screwing the sensor into the process. O-rings are consumable parts which must be exchanged regularly (at least once per year). Even when all required safety measures have been complied with, potential risks still exist with respect to leaks or mechanical damage to the armature. Wherever there are seals or screws, gases or liquids can leak out undetected.

Before removing the sensor from its measurement setup, always make sure that no process medium can be accidentally spilled.

The built-in temperature sensor can only be used for monitoring the sensor conditions, not for controlling the process temperature.

4. Initial operation

This CONDUCELL sensor has been carefully tested and is ready for use. Always check a sensor for possible defects after first unpacking it. In the unlikely event that you find a damaged sensor, return immediately the CONDUCELL in its original packing to your HAMILTON representative.

1: Socket head with VP8 connector, 2: ARC logo, serial and part number, 3: PG 13.5 thread, 4: O-ring, 5: Sensor shaft 6: Sensor tip with current and potential electrodes.

The CONDUCELL 4USF ARC sensors are delivered directly from the factory already pre-calibrated (according to the operational specifications). The integrated 4–20 mA analog interface and RS485 digital interface (Modbus RTU) are configured according to factory defaults. You can find full details, including serial number and most important specifications, on the certificate provided with each sensor.

ATTENTION!


To avoid contamination problems, make sure the sensor is always stored with a protection cap on the sensor tip.

To be sure to avoid electrical damage to the sensor, carefully follow all the instructions in the section entitled 'Electrical Connection'.

Before application of the sensor for measurement, monitoring or regulation, be sure to first check its configuration by means of an operational test.

5. Electrical connection

The CONDUCELL 4USF ARC sensor is fitted with a VP8 socket head. The eight golden contacts are denoted as pin A... pin H. For easy identification of each pin the head has a mark between pin A and pin B.

For the easiest and safest connection of CONDUCELL 4USF ARC sensors, always use HAMILTON VP8 cables. These are available in a range of different lengths.

VP pin	Function
A	4–20 mA interface (mA interface #2).
B	4–20 mA interface (mA interface #1).
C	Power supply: +24 VDC (7 to 30 VDC).
D	Power supply: Ground.
G	RS485 (A).
H	RS485 (B).

Electrical connection of the 4–20 mA current interfaces

The 4–20 mA interface enables direct connection of the CONDUCELL 4USF ARC sensor to a data recorder, indicator, control unit or PLC with analog I/O. Apart from the two wires used for 4–20 mA current loop no additional equipment is required for analog signaling.

CONDUCELL 4USF ARC sensors are delivered from the factory with the analog 4–20 mA interface active (see the included sensor certificate).

When using the 4–20 mA interface, pins have the following designations with respect to VP cable conductor colors:

CONDUCELL 4USF ARC	VP pin	VP8 cable
4–20 mA three-wire interface, functions as a current sink. It regulates the input current according to the sensor measurements. This interface needs a separate power supply (pin B or C). Factory default is assigned to a temperature measurement.	A	Coaxial core black transparent.
4–20 mA two-wire interface, functions as a current sink. It regulates the input current according to the sensor measurements. This interface can be powered directly from 4–20 mA two-wire current loop. Nominal power of 60 mW must be provided. Factory default is assigned to a conductivity measurement.	B	Coaxial shield black.
Power supply: +24 VDC (7 to 30 VDC). Max. power consumption: 150 mW.	C	Coaxial core red transparent.
Power supply: Ground.	D	Coaxial shield red.

The 4–20 mA interface is configured at the factory with the range, value and unit for the measurements as indicated on the certificate. Follow the instructions in the section entitled 'Configuration and monitoring of the sensor' to adjust the sensor according to the requirements of your application.

Examples of circuit arrangement


Fig. A: Two-wire loop wiring diagram for the 4–20 mA interface (mA interface #1). In this wiring scheme the power is not supplied to the sensor VP pin C, wherefore the wiring is not applicable to sensor with ARC Wi Sensor Adapter and VISICAL.


Fig. B: Three-wire loop wiring diagram for the 4–20 mA interfaces. The figure represents both 4–20 mA interfaces at pin A and pin B.


Fig. C: The safest form of wiring, using an isolation amplifier. The figure represents both 4–20 mA interfaces at pin A and pin B. (For detailed technical advice, please contact the technical support at HAMILTON.)

Electrical connection for the digital RS485 interface

The digital RS485 interface enables communication with CONDUCELL 4USF ARC sensor for performing measurements, for calibrating the sensor and for changing the sensor's configuration parameters. CONDUCELL 4USF ARC sensors are always connected to digital controlling devices as a Modbus slave. To function, they require a power supply (VP8 pins C and D, see below). The section entitled 'Configuration and monitoring of the sensor' describes operation in digital mode.

By using the correct access password, the ARC System operator can adapt ARC sensors to many tasks by:

- Selecting the 4–20 mA interface.
- Scaling (configuring) the 4–20 mA interface.
- Selecting the measured parameter:
 - Specific conductivity: $\mu\text{S}/\text{cm}$; mS/cm .
 - Temperature T: $^{\circ}\text{C}$; K.
- Define the Temperature compensation coefficient.

In addition, operators can read sensor information from the RS485 interface such as:

- The sensor's serial number, part number (P/N) and manufacturing number (WO).
- The sensor's firmware version.
- The sensor's status (e.g. operation hours, number of cleanings and sterilizations, warnings and errors).

Additional information:

The Modbus RTU communication protocol corresponds to the Modbus-IDA standard (see <http://www.modbus.org>). CONDUCELL 4USF ARC uses an open register set developed by HAMILTON. Additional information about the register content and structure can be found at <http://www.hamiltoncompany.com>.

ATTENTION!

Because all sensors are delivered with factory-default settings, each sensor must be configured for its specific application before first use (see the section entitled 'Configuration of the sensor').

The pins for digital the RS485 interface have the following designation with respect to VP cable conductor colors:

CONDUCELL 4USF ARC	VP pin	VP8 cable
Power supply: +24 VDC (7 to 30 VDC) Max. power consumption 150 mW.	C	Coaxial core red transparent.
Power supply: Ground.	D	Coaxial shield red.
RS485 (A).	G	Yellow.
RS485 (B).	H	Brown.

In an electromagnetically noisy environment, it is advisable to connect the VP cable shield to the ground. This significantly improves noise immunity and signal quality.

Example of a circuit arrangement


Fig. D: Wiring diagram for the RS485 interface.


Fig. E: Multi-drop bus wiring for the Modbus two-wire mode. Each sensor functions as a Modbus slave.

NOTE!

In the connection scheme shown above, each sensor must have the unique Modbus device address for proper communication.

The serial Modbus connection between the RS485 port of the master and the corresponding interfaces of the sensors has to be ensured according to the EIA/TIA RS485 standard. Only one sensor can communicate with the master at any time.

6. Configuration and monitoring of the sensor

Two options are available for the configuration and monitoring of a sensor:

1. Personal computer or notebook. The following additional equipment and software is also required:
 - HAMILTON USB-RS485 Modbus converter (P/N: 242411)
 - HAMILTON freeware 'ARC Sensor Configurator' available at <http://www.hamiltoncompany.com> (follow the instructions in the Configurator's user guide for the installation and operation).
 - Demo Cable (P/N: 355194). This cable includes a power adapter to supply the sensor with operation power, and a plug to connect the two RS485 conducting wires (yellow and brown) to USB-RS485. If you use a standard VP8 cable, you must supply the sensor with an external power source (pin C: 24 VDC; pin D: Ground).
2. HAMILTON ARC View Handheld Package (P/N: 242180). The handheld package represents an ideal desktop solution for the ARC sensor management. The ARC View Handheld included in the package is a compact mobile wireless device with long battery life and broad functionality. When using the Handheld as a mobile wireless device, each ARC sensor requires an ARC Wi Sensor Adapter (P/N: 242170).

The digital RS485 interface of the CONDUCELL 4USF ARC sensor can be accessed by operators working within a three-tier hierarchy. Operator hierarchy levels and factory default passwords are shown in the table below.

Operator status	Operator Level	Password
User	U	Not required
Administrator	A	18111978
Specialist	S	16021966

The Users can read basic data from sensor. Administrators can also calibrate sensors. Specialists can calibrate and configure sensors, and can see all data.

7. Preparation for the measurement

The cell constant and the linearity of the sensor can vary with the fitting situation.

A symmetrical setup is recommended in order to ensure a high degree of linearity. The diameter of the measuring equipment should be min. 30mm. and the distance from the front side - min. 10 mm. To achieve high precision, the cell constant should be calibrated in the final setup.

Prepare the sensor for measurements as follows:

1. Remove the protective caps from the VP head and from the sensor tip.
2. Make sure that the sensor is configured as required. If in doubt, test as described in the section entitled 'Electrical connection: RS485 Modbus interface'.
3. Installation of the sensor (thread: PG 13.5).

4. Attach the CONDUCELL 4USF ARC sensor according to the section 'Electrical Connection' in the desired configuration (analog 4–20 mA interface, digital RS485 interface or both).

The sensor signal stabilizes itself within few minutes. The sensor is factory pre-calibrated and ready for use. However, if an on-site calibration is necessary, a calibration can be performed via RS485 by means of the ARC View Handheld, personal computer or VISICAL (P/N: 242410).

8. Removal of the sensor

Before removing any sensor from a process, make sure that no process medium can escape during removal.

Unscrew the PG 13.5 thread connection and pull out the sensor.

9. Sterilization, autoclaving, CIP procedures

CONDUCELL 4USF ARC sensors are carefully designed to withstand the cleaning procedures normally used in biotechnological and general chemistry industries without the need for special precautions. Nevertheless, as is the case with all sensors, it is seen that frequent cleanings or sterilizations with steam or hot and aggressive solutions lead to shortened life span of the sensor. Although for a number of reasons (variations in composition of cleaning agent, temperature gradient and so on) it is not possible to give general life span figures. With steam sterilizations and autoclaving for 30 min at a temperature of 125°C, a life span of more than 50 cycles per exchangeable sensor is common.

The contacts must be clean and dry before sensor is connected to the cable.

ATTENTION!

The sensor can measure and communicate over digital RS485 interface up to a temperature of 130°C. However, the 3.5 mA signal is set on both 4–20 mA interfaces as the process temperature exceeds 110°C. The analog measurements will be continued when the temperature goes back below 110°C.

NOTE!

The cell constant of the sensor can change by corrosion of the electrodes in aggressive media, high temperature or by contamination of the electrodes during the sensors lifetime. A quality indicator of the CONDUCELL 4USF ARC sensor shows deviation of the cell constant. The quality indicator status is updated automatically after each calibration.

10. Testing and maintenance

Calibration

The CONDUCELL 4USF ARC sensors provide two kinds of sensor calibration: automatic standard calibration, and product calibration.

The automatic standard calibration is performed over RS485 interface using the ARC View Handheld, VISICAL or the ARC Sensor Configurator freeware (see the sections entitled 'Configuration and monitoring of the sensor' and 'Accessories'). For the product calibration use the ARC View Handheld, or the ARC Sensor Configurator freeware (VISICAL is not supporting the product calibration.)

Automatic standard calibration

CONDUCELL 4USF ARC sensor utilizes a one-point calibration procedure with automatic standard recognition. The conductivity values and temperature dependence data for a number of predefined conductivity standards are stored in the sensor. Therefore the sensor recognizes a standard, examines the correctness and stability of the conductivity and temperature signals, and performs the calibration.

Actions:

1. Set an appropriate operator level (Administrator or Specialist).
2. Check that the correct set of the calibration standards is selected.
3. Immerse sensor into corresponding conductivity standard.
4. Let the system equilibrate. Guarantee stable conditions for at least three minutes.
5. Execute the calibration.
6. If all conditions are met, the sensor confirms the calibration immediately.

Three sets of the conductivity calibration standards are provided with each CONDUCELL 4USF ARC sensor:

- HAMILTON
- REAGECON
- KCI solutions

However, the sensors are pre-calibrated at the factory with HAMILTON 1413 $\mu\text{S}/\text{cm}$ conductivity standard at room temperature. The following standards from HAMILTON standards set are defined as factory default for calibration with automatic standard recognition: 1.3 $\mu\text{S}/\text{cm}$, 5 $\mu\text{S}/\text{cm}$, 15 $\mu\text{S}/\text{cm}$, 84 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$, 12880 $\mu\text{S}/\text{cm}$, and 100 mS/cm .

The use of standards with other conductivities will yield an error and reject of the calibration.

NOTE!

To ensure optimal measurement accuracy always use conductivity standard with value which fits most to the measuring range of your application.

HINT: The concept behind CONDUCELL 4USF ARC sensors enables calibration and configuration in the lab before use in the process control. Another calibration for the installation in the process setup is not required. However, if the propagation of the electrical field is influenced by the mounting situation an adjustment of the cell constant or a calibration on-site might be necessary.

NOTE!

The ARC View Handheld supports automatic calibration for a predefined set of conductivity standards. Use ARC Sensor Configurator software to perform manual calibration or to configure the set of standards.

Product calibration.

The product calibration is an in-process calibration procedure in order to adjust the measurement to specific process conditions, or in case the sensor cannot be removed for the standard calibration.

Product calibration is an additional calibration procedure to an automatic standard calibration. Product calibration adapts the cell constant to the process conditions in force at the time of product calibration. In order to restore the original cell constant, the product calibration can be cancelled at any time. A new standard calibration cancels a product calibration as well.

Actions:

1. Set an appropriate operator level (Administrator or Specialist).
2. Perform an initial measurement while taking a sample from the process. The data of the initial measurement are stored in the sensor.
3. Perform a laboratory measurement of the sample at the same temperature as it was measured in the process. The laboratory reference device must have the same settings as the CONDUCELL 4USF ARC sensor regarding temperature compensation and reference temperature.
4. Assign the laboratory value to the value of the initial measurement. The product calibration is accepted and instantaneous active if the corresponding cell constant is changed not more than 70% of its original value.

NOTE!

The reference temperature for the temperature compensation function of CONDUCELL 4USF ARC sensor is 25 °C..

Self-diagnosis functions

CONDUCELL 4USF ARC sensors provide a self-diagnosis functionality to detect and identify the most common sensor malfunctions. Both interfaces can be used for warning and error messaging. The analog 4–20 mA interface can be configured according to the NAMUR recommendations to indicate an abnormal event. The RS485 interface provides a variety of the indications based on the fault code.

Use the ARC View Handheld and ARC Wi Sensor Adapter for monitoring of the sensor status and troubleshooting. These devices detect and highlight the sensor condition according to the fault code. The ARC View Handheld displays the text messages that are appropriate for the type of malfunction.

The following types of messages are provided by the self-diagnosis function:

- Warning (calibration):
 - Calibration recommended.
 - Calibration upper/lower, out of range, or unstable.
- Error (failure):
 - Conductivity sensor defect (reading failure).
 - Temperature upper/lower, out of range.
 - Temperature sensor defect (reading failure).
 - Electrodes resistance out of range (for example in air).
 - Sensor quality low (calibration is not possible).
 - Internal communication failure.

11. Disposal

The design of HAMILTON sensors optimally considers environmental compatibility. In accordance with the EC guideline 2002/96/EG HAMILTON sensors that are worn out or no longer required must be sent to a dedicated collection point for electrical and electronic devices, alternatively, must be sent to HAMILTON for disposal. Sensors must not be sent to an unsorted waste disposal point.

12. Parts and Accessories

ARC View Handheld with Dock.

ARC Wi Sensor Adapter for wireless communication.

VISICAL.

USB-RS485 Modbus Converter.

Sensor Cable VP8.

Order number	Description
242410	VISICAL
242411	USB-RS485 Modbus converter
242180	ARC View Handheld Package
355194	Demo Cable (1m, open end, with plug power pack)
242170	ARC Wi Sensor Adapter
238999-2394	Sensor cable VP8, 1m
238999-1953	Sensor cable VP8, 3m
238999-2395	Sensor cable VP8, 5m
238999-2396	Sensor cable VP8, 10m
238999-2403	Sensor cable VP8, 15m
238999-2505	Sensor cable VP8, 20m

13. Technical data

- CONDUCELL 4USF sensor with integrated electronics.
Functionality includes measurement and self-diagnosis.
- Steam sterilizable, autoclavable and compatible with Cleaning In Place (CIP).
- Measurement principle: 4 electrode contacting measurement.
- Shaft has a diameter of 12 mm and a PG 13.5 thread.
- Different lengths of shaft available, from 120 mm to 425 mm.
- Operating temperature range: -20 to 130°C.
- Process pressure range: -1 to 20 bar / 290 psi (relative).
- Range of measurement: 1 to 500'000 µS/cm.
- Measured values can be configured by software according to:
 - Conductivity: µS/cm, mS/cm.
 - Temperature °C, K.
- Medium affected materials: see certificate.
- Electrical connection: VP8 socket head.
- 7 to 30 VDC operating voltage; max. 150 mW.
- Two scalable 4–20 mA current interfaces (two-wire) for conductivity and temperature signals.
- Digital RS485 Interface (two-wire):
 - Protocol: Modbus RTU; maximal 31 addresses.
 - Baud rate: 4800, 9600, 19200, 38400, 57600, 115200 bd.

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Copyright © 2010 HAMILTON Bonaduz AG, Bonaduz Schweiz. Alle Rechte vorbehalten. Die Reproduktion irgendeines Teils dieses Dokuments in jeder beliebigen Form ist ohne die ausdrückliche schriftliche Zustimmung der HAMILTON Bonaduz AG untersagt.