

## HAMILTON GASTIGHT® INSTRUMENT SYRINGES

Congratulations! You have purchased the finest quality precision syringe available today. We at Hamilton Company combine top quality materials with skilled workmanship, ensuring the highest possible performance level of every precision fluid device we manufacture. With proper care and handling, GASTIGHT Instrument Syringes will provide unsurpassed performance in precision liquid handling year after year.

Hamilton GASTIGHT Instrument Syringes are precision, measuring instruments. For prolonged syringe life and to obtain the maximum benefits of use, a few helpful tips should be observed:

### ASSURED ACCURACY AND PRECISION

- When initially inserting a new plunger into a GASTIGHT syringe barrel, wet (lubricate) the tip with deionized water or another solvent compatible with system fluid.
- For extended plunger tip life, minimize the use of a "dry" syringe.
- Every GASTIGHT Instrument Syringe is handcrafted to ensure only the highest precision liquid transfers. Barrels and plungers within a given volume range are generally interchangeable and may be purchased separately for field repair and replacement.
- If the plunger is accidentally withdrawn completely from the syringe barrel, wipe it carefully with a lint-free tissue and re-wet it before reinserting into the barrel. Be careful of the plunger tip since any physical abrasions, scratches or oil from one's fingers may cause the plunger to leak once re-assembled.

### CLEANING

- The life of your GASTIGHT syringe is directly related to its cleanliness!
- To clean syringes, it is best to use solvents known to be effective in solvating the system fluid.
- High quality grade water and acetone prove to be good rinses.

- It is recommended that the plunger not be removed from the barrel of an instrument syringe for cleaning. Precipitants can form a ring of contamination in the barrel at the end of the fill stroke. Moving the plunger past this ring can damage the tip.

### PRESSURE PARAMETERS

Syringe Model	Pressure (psi)
1702-1710	1000
1725.1-1750	500
1001-1010	200
1025	100

### TEMPERATURE PARAMETERS

- Operating: 10° C to 40° C.
- Storage: -20° C to 70° C.
- Even though only the highest quality glass is used for GASTIGHT syringe barrels, rapid changes in temperature should be avoided.

### WARRANTY STATEMENT

Hamilton Company unconditionally guarantees its products to be free of defects in materials and workmanship. Any product which fails due to such defects will be repaired or replaced at our discretion without cost, provided the device is returned with an explanation. It is the responsibility of the purchaser to determine the suitability of application and material compatibility of the products based on the published specifications of the products.



**HAMILTON**  
THE MEASURE OF EXCELLENCE™

[www.hamiltoncompany.com](http://www.hamiltoncompany.com)

**Hamilton Company**

4970 Energy Way  
Reno, Nevada 89502 USA  
Toll-Free: 800-648-5950  
Telephone: +1-775-858-3000  
Fax: +1-775-856-7259  
e-mail: [sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

**Hamilton Bonaduz AG**

Via Crusch 8  
CH-7402  
Bonaduz/Switzerland  
Telephone: +41-(0)81-660-60-60  
Fax: +41-(0)81-660-60-70  
e-mail: [marketing@hamilton.ch](mailto:marketing@hamilton.ch)

**Hamilton GmbH**

Daimlerweg 5A  
64293 Darmstadt/Germany  
Telephone: +49-(0)6151-66706-0  
Fax: +49-(0)6151-66706-66  
e-mail: [sales@hamiltongmbh.de](mailto:sales@hamiltongmbh.de)

**Hamilton Northern Europe**

Unit 2, Enterprise Way  
Aston Science Park  
Birmingham, B7 48H, U.K.  
Telephone: +44-(0)121 2600301  
Fax: +44-(0)121 2600302  
e-mail: [sales@hamiltongb.co.uk](mailto:sales@hamiltongb.co.uk)

**Hamilton France S.A.R.L.**

Parc de Moulin de Massy  
37 rue du Saule Trapu  
91300 MASSY/France  
Telephone: +33-(0)1-69-75-16-16  
Fax: +33-(0)1-60-11-57-16  
e-mail: [Hamilton.France@wanadoo.fr](mailto:Hamilton.France@wanadoo.fr)

**Quality Hamilton Products:**

MICROLITER™ Syringes  
GASTIGHT® Syringes  
Chromatography Syringes  
Syringes for Life Science  
Instrument Syringes  
SoftGrip™ Pipettes  
SoftTouch™ Electronic Pipettes  
SoftAide™ Pipette Controllers  
Miniature Valves  
Modular Valve Positioner  
Laboratory Fittings, Adapters & Tubing  
Precision Syringe Pumps  
MICROLAB® Diluters & Dispensers  
SoftTop™ Quik Dispensers  
Electrochemical Sensors  
DURACAL™ Buffer Solutions  
HPLC Columns & Resins

**Laboratory Automation for:**

Drug Discovery  
Genomics  
Proteomics  
Forensic  
In Vitro Diagnostics

## RETURN OF GOODS

Hamilton Company's return and repair policy is written to protect its employees from potentially hazardous materials (e.g. serum, radioactive materials, carcinogenic chemicals, etc.) or any substance that may cause them partial or permanent disability during the inspection or repair process. In returning product, the customer acknowledges that the product is free from any hazardous materials. Furthermore, the customer assumes responsibility if the returned product is determined to be hazardous.

**Syringes and needles manufactured by Hamilton Company are intended for scientific research and laboratory use only and are not intended for human in vivo use.**

Document No. 53036-01 REV. B

© Copyright Hamilton Company 4/06

Printed in the U.S.A.