# VISIFERM™ DO SENSORS

Operating Instructions

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### **Important note**

Copyright © 2010 HAMILTON Bonaduz AG, Bonaduz Switzerland. All rights reserved. The reproduction of any part of this document in any form is forbidden without the express written agreement of HAMILTON Bonaduz AG.

Contents of this manual can be modified without previous announcement. Technical modifications reserved. Greatest possible care was used on the correctness of the information in this manual. If errors should be discovered nevertheless, HAMILTON Bonaduz AG is pleased to be informed about it. Regardless of this, HAMILTON Bonaduz AG cannot assume liability for any errors in this manual or for their consequences.

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<b>1. INTRODUCTION</b>	<b>4</b>
<b>2. INTENDED USE</b>	<b>4</b>
<b>3. SAFETY INSTRUCTIONS</b>	<b>5</b>
<b>4. INITIAL OPERATION</b>	<b>5</b>
<b>5. ELECTRICAL CONNECTION</b>	<b>6</b>
ELECTRICAL CONNECTION OF THE 4–20 mA CURRENT INTERFACES	6
EXAMPLES OF CIRCUIT ARRANGEMENT	7
ELECTRICAL CONNECTION FOR THE ECS MODE	8
EXAMPLE OF CIRCUIT ARRANGEMENT	9
ELECTRICAL CONNECTION FOR THE DIGITAL RS485 INTERFACE	9
EXAMPLE OF CIRCUIT ARRANGEMENT	10
<b>6. CONFIGURATION AND MONITORING OF THE SENSOR</b>	<b>11</b>
<b>7. PREPARATION FOR THE MEASUREMENT</b>	<b>12</b>
<b>8. REMOVAL OF THE SENSOR</b>	<b>12</b>
<b>9. STERILIZATION, AUTOCLAVING, CIP PROCEDURES</b>	<b>13</b>
<b>10. TESTING AND MAINTENANCE</b>	<b>13</b>
CALIBRATION	13
CALIBRATION IN ECS MODE	15
CHANGING THE SENSOR CAP	15
SELF-DIAGNOSIS FUNCTIONS	16
<b>11. DISPOSAL</b>	<b>16</b>
<b>12. PARTS AND ACCESSORIES</b>	<b>17</b>
<b>13. TECHNICAL DATA</b>	<b>18</b>
<b>BEDIENUNGSANLEITUNG DEUTSCH</b>	<b>20</b>

# Operating Instructions VISIFERM™ DO

## 1. Introduction

This manual refers to the VISIFERM™ DO range of optical, sterilizable oxygen sensors from HAMILTON Bonaduz AG.

Designation	Order number P/N
VISIFERM™ DO 120 (4-20 mA mode)	242 450-01
VISIFERM™ DO 120 (ECS mode)	242 450-02
VISIFERM™ DO 160 (4-20 mA mode)	242 451-01
VISIFERM™ DO 160 (ECS mode)	242 451-02
VISIFERM™ DO 225* (4-20 mA mode)	242 452-01
VISIFERM™ DO 225* (ECS mode)	242 452-02
VISIFERM™ DO 325 (4-20 mA mode)	242 453-01
VISIFERM™ DO 325 (ECS mode)	242 453-02
VISIFERM™ DO 425 (4-20 mA mode)	242 454-01
VISIFERM™ DO 425 (ECS mode)	242 454-02

\*The VISIFERM™ DO 225 has, in reality, a shaft length of 215 mm. This is to ensure optimal rinsing characteristics in replaceable armatures such as RETRACTEX.

Hamilton VISIFERM™ DO sensors are quality products manufactured according to the very latest research findings. Follow the instructions given here to be sure of optimal safety and durability.

**IMPORTANT:** These instructions must be read, understood and followed by all staff using VISIFERM™ DO sensors. HAMILTON can assume no responsibility for damage and operational disruption arising from failure to observe these instructions.

## 2. Intended use

VISIFERM™ DO sensors were developed for the measurement of dissolved oxygen partial pressure, together with the following derived measurement parameters:

- Percent by volume oxygen.
- Percent air saturation by oxygen.
- Concentration of oxygen in liquids.

The main application for VISIFERM™ DO sensors is the in-line measurement under demanding conditions in biotechnology. The main characteristic that makes VISIFERM sensors ideal for this application is their remarkable long-term stability, even after frequent sterilizations and autoclavings.

Nevertheless, in many other applications such as those found in the chemical industry, air surveillance, fish farming, water management, and sewage management, VISIFERM™ DO sensors have also already proven their value.

During development, special attention was paid to an optimum sanitary

design. All materials in contact with the solution meet FDA requirements.

VISIFERM™ DO sensors provide electrical interfaces: ECS (classical electrochemical sensor), standard analog (4–20 mA) and digital RS485 (Modbus RTU). These are built into each sensor and are supported directly from the sensor head. The analog 4–20 mA and digital RS485 interfaces do not require any additional equipment such as amplifier or transmitter.

Beside different interfaces, the VISIFERM™ DO sensor also contains a temperature sensor (NTC 22 KOhm). This temperature sensor only has to be used for the compensation of temperature of the oxygen signal in ECS mode (simulation of a classical amperometric oxygen sensor by means of VISIFERM™ DO). The temperature sensor must not be used for the regulation of the process temperature.

### 3. Safety instructions

VISIFERM™ DO sensors must be used for their intended applications, and in optimum safety and operational conditions. The specifications (such as temperature or pressure) defined in the section entitled 'Technical data' must not be exceeded under any circumstances. Potential hazards can exist if the sensor is not operated correctly or appropriately.

Assembly and maintenance must be performed only by trained personnel.

Make sure that the PG 13.5 thread and the O-ring are not damaged when screwing the sensor into the process. O-rings are consumable parts which must be exchanged regularly (at least once per year). Even when all required safety measures have been complied with, potential risks still exist with respect to leaks or mechanical damage to the armature. Wherever there are seals or screws, gases or liquids can leak out undetected.

Before removing the sensor from its measurement setup, always make sure that no process medium can be accidentally spilled.

The built-in temperature sensor can only be used for monitoring the sensor conditions, not for controlling the process temperature.

### 4. Initial operation

This VISIFERM™ DO sensor has been carefully tested and is ready for use. Always check a sensor for possible defects after first unpacking it. In the unlikely event that you find a damaged sensor, return immediately the VISIFERM™ DO in its original packing to your HAMILTON representative.

1: Socket head with VP8 connector, 2: Key areas with serial number, 3: PG 13.5 thread, 4: O-ring, 5: Sensor shaft with stainless steel heat number, 6, 7: Sensor cap with oxygen-sensitive silicone luminophore (sensory element), sensor cap serial number, and stainless steel heat number.

The VISIFERM™ DO sensors are delivered directly from the factory already pre-calibrated (according to the operational specifications). The integrated analog ECS and 4–20 mA interfaces and RS485 digital interface (Modbus RTU) are configured according to factory defaults. You can find full details, including serial number and most important specifications, on the certificate provided with each sensor.

**ATTENTION !**

*In the ECS mode the sensor can be demolished if you try to operate the sensor as a 4-20 mA sensor!*

**ATTENTION!**

*To avoid humidity problems, make sure that the sensor cap is always attached firmly to the sensor shaft, and that the O-ring between the shaft and cap is undamaged.*

To be sure to avoid electrical damage to the sensor, carefully follow all the instructions in the section entitled 'Electrical Connection'.

Before application of the sensor for measurement, monitoring or regulation, be sure to first check its configuration by means of an operational test.

**5. Electrical connection**

The VISIFERM™DO sensor is fitted with a VP8 socket head. The eight golden contacts are denoted as pin A... pin H. For easy identification of each pin the head has a mark between pin A and pin B.

For the easiest and safest connection of VISIFERM™DO sensors, always use HAMILTON VP8 cables. These are available in a range of different lengths.

VP pin	Function
A	cathode (only in ECS mode)
B	- 4 – 20 mA interface (mA interface #1). - ECS mode: anode. Never connect in the ECS mode with +2 V or more!
C	Power supply: +24 VDC (7 to 30 VDC). Start-up power: 1 W. Continuous power consumption: 600 mW.
D	Power supply: Ground.
E	Temperature sensor NTC 22k for ECS mode
F	Temperature sensor NTC 22k for ECS mode
G	RS485 (A).
H	RS485 (B).
Casing	Sensor shaft shielding. The outside shielding of the VP cable should be connected to power supply ground (Pin D).

**Electrical connection of the 4–20 mA current interfaces**

The 4–20 mA interface enables direct connection of the VISIFERM™DO sensor to a data recorder, indicator, control unit or PLC with analog I/O. Apart from the two wires required for the (temperature compensated) 4–20 mA oxygen signal only a power supply is required for the VISIFERM™ DO sensor (VP8 pins C and D, see below).

**! ATTENTION !**

*Ensure that the sensor is adjusted at the 4-20 mA mode! If the sensor is configured in the ECS mode, and if you use it as a 4-20 mA sensor, then this operation can result in irreparable damages in the sensor!*

When using the 4–20 mA interface, pins have the following designations with respect to VP cable conductor colors:

VISIFERM™ DO	VP pin	VP8 cable
4–20 mA three-wire interface, functions as a current sink. It regulates the input current according to the sensor measurements. This interface needs a separate power supply (pin C). Factory default is assigned to an oxygen measurement.	B	Coaxial shield black.
Power supply: +24 VDC (7 to 30 VDC). Start-up power: 1 W. Continuous power consumption: 600 mW.	C	Coaxial core Red transparent.
Power supply: Ground.	D	Coaxial shield red.
Sensor shaft (connected to power supply ground).	Shield	Cable shield green-yellow.

In an electromagnetically noisy environment, it is advisable to assign the sensor's shaft and/or VP cable shield to ground or earth. This significantly improves noise immunity and signal quality.

The 4–20 mA interface is configured at the factory with the range, value and unit for the measurements as indicated on the certificate. Follow the instructions in the section entitled 'Configuration and monitoring of the sensor' to adjust the sensor according to the requirements of your application.

**Examples of circuit arrangement**


Fig. A: Three-wire loop wiring diagram for the 4–20mA interface.


Fig. B: The safest form of wiring, using an isolation amplifier. The figure represents 4–20 mA interface. (For detailed technical advice, please contact the technical support at HAMILTON.)

### Electrical connection for the ECS mode

The ECS mode enables the simulation of an electrochemical sensor. Thus a VISIFERM™ DO sensor can be connected to classical measuring devices instead of amperometric oxygen sensors (Clark cells). Furthermore only the power supply of the VISIFERM™ DO sensor is necessary.

#### ATTENTION !

*Do not apply any high voltage (max. 2 VDC) at pin B (anode)! This can result in a destruction of the sensor in ECS mode! Note: Only in 4-20 mA mode a high voltage (max. 24 VDC) may be applied in order to operate the current interface!*

When using the ECS interface, pins have the following designations with respect to VP cable conductor colors:

VISIFERM™ DO	VP pin	VP8 cable
cathode	A	Coaxial core black transparent
anode <b>ATTENTION:</b> In ECS mode never connect with a potential higher as +2 V!	B	Coaxial shield black
power supply: +24 VDC (7 to 30 VDC), max. 1000 mW	C	Coaxial core red transparent
power supply: ground	D	Coaxial shield red
NTC 22 KOhm	E	White
NTC 22 KOhm	F	Green
sensor shaft (connect with the mass of the power supply)	shield	Cable shield green-yellow

In an electromagnetically noisy environment, it is advisable to assign the sensor's shaft and/or VP cable shield to ground or earth. This significantly improves noise immunity and signal quality.

The NTC temperature sensor attached to the pins E and F is isolated from the remainder of the integrated electronics and is used for the temperature compensation of the oxygen signal in the measuring device. If the measuring

device is not designed for typical sterilizable oxygen sensors, then the measuring device probably expects another type of temperature sensor. Therefore the measuring device neither can compensate the temperature depending oxygen signal nor display the temperature correctly.

Usually classical sensors are operated with a polarization voltage between anode and cathode. This polarization voltage is supplied by the measuring device. VISIFERM™ DO can be operated with polarization voltages usual for electrochemical sensors. The sensor is optimized for a polarization voltage of -675mV.

For adjustment to different measuring devices and/or for simulation of different amperometric sensors the current can be adjusted between 0 and 500 nA.

### Example of circuit arrangement


Fig. C: Wiring diagram for the ECS interface.

### Electrical connection for the digital RS485 interface

The digital RS485 interface enables communication with VISIFERM™ DO sensor for performing measurements, for calibrating the sensor and for changing the sensor's configuration parameters. VISIFERM™ DO sensors are always connected to digital controlling devices as a Modbus slave. To function, they require a power supply (VP8 pins C and D, see below). The section entitled 'Configuration and monitoring of the sensor' describes operation in digital mode.

By using the correct access password, the operator can adapt sensor to many tasks by:

- Switching between the 4–20 mA and ECS interface mode.
- Scaling (configuring) the 4–20 mA, or ECS interfaces.
- Selecting the measured parameter:
  - DO: % oxygen air saturation (%-sat), volume % of oxygen (%-vol), mg/l or ppm, ug/l or ppb.
  - Temperature T: °C.
- Compensating for the process pressure, for calculation of volume percentage values (%-sat, %-vol).
- Compensating for the salinity of the measured solution, for calculation of the concentration values (mg/l or ppm, ug/l or ppb).

In addition, operators can read sensor information from the RS485 interface such as:

- The sensor's serial number, part number (P/N) and manufacturing number (WO).
- The sensor's firmware version.
- The sensor's status (e.g. operation hours, number of cleanings and sterilizations, warnings and errors).

Additional information:

The Modbus RTU communication protocol corresponds to the Modbus-IDA standard (see <http://www.modbus.org>). VISIFERM™ DO uses an open register set developed by HAMILTON. Additional information about the register content and structure can be found at <http://www.hamiltoncompany.com>.

**ATTENTION!**

*Because all sensors are delivered with factory-default settings, each sensor must be configured for its specific application before first use (see the section entitled 'Configuration and monitoring of the sensor').*

The pins for digital the RS485 interface have the following designation with respect to VP cable conductor colors:

VISIFERM™ DO	VP pin	VP8 cable
Power supply: +24 VDC (7 to 30 VDC), max. power consumption 1 W.	C	Coaxial core red transparent.
Power supply: Ground.	D	Coaxial shield red.
RS485 (A).	G	Yellow.
RS485 (B).	H	Brown.
Sensor shaft (connected to the ground of the power supply).	Shield	Cable shield green-yellow.

In an electromagnetically noisy environment, it is advisable to connect the VP cable shield to the ground. This significantly improves noise immunity and signal quality.

**Example of circuit arrangement**


Fig. D: Wiring diagram for the RS485 interface.


Fig. E: Multi-drop bus wiring for the Modbus two-wire mode. Each sensor functions as a Modbus slave.

**NOTE!**

*In the connection scheme shown above, each sensor must have the unique Modbus device address for proper communication.*

The serial Modbus connection between the RS485 port of the master and the corresponding interfaces of the sensors has to be ensured according to the EIA/TIA RS485 standard. Only one sensor can communicate with the master at any time.

## 6. Configuration and monitoring of the sensor

Two options are available for the digital configuration and monitoring of a sensor:

1. Personal computer or notebook. The following additional equipment and software is also required:
  - HAMILTON USB-RS485 Modbus converter (P/N: 242411)
  - HAMILTON freeware 'ARC Sensor Configurator' available at <http://www.hamiltoncompany.com> (follow the instructions in the Configurator's user guide for the installation and operation).
  - Demo Cable (P/N: 355194). This cable includes a power adapter to supply the sensor with operation power, and a plug to connect the two RS485 conducting wires (yellow and brown) to USB-RS485. If you use a standard VP8 cable, you must supply the sensor with an external power source (pin C: 24 VDC; pin D: Ground).
2. HAMILTON ARC View Handheld Package (P/N: 242180). The handheld package represents an ideal desktop solution for the ARC sensor management. The ARC View Handheld included in the package is a compact mobile wireless device with long battery life and broad functionality. When using the Handheld as a mobile wireless device, each ARC sensor requires an ARC Wi Sensor Adapter (P/N: 242170).

The digital RS485 interface of the VISIFERM™ DO sensor can be accessed by operators working within a three-tier hierarchy. Operator hierarchy levels and factory default passwords are shown in the table below.

Operator status	Operator Level	Password
User	U	Not required
Administrator	A	18111978
Specialist	S	16021966

The Users can read basic data from sensor. Administrators can also calibrate sensors. Specialists can calibrate and configure sensors, and can see all data.

## 7. Preparation for the measurement

In contrast to electrochemical oxygen sensors, VISIFERM™ DO sensors work independently of installation position orientation.

Prepare the sensor for measurements as follows:

1. Remove the protective caps from the VP head and from the sensor cap.
2. Make sure that the sensor is configured as required. If in doubt, test as described in the section entitled 'Electrical connection: RS485 Modbus interface'.
3. Check the condition of the sensor and sensor cap.
4. Make sure that the sensor cap is attached firmly to the sensor shaft, and that the O-ring between the shaft and cap is undamaged.
5. Installation of the sensor (thread: PG 13.5).
6. Attach the VISIFERM™ DO sensor according to the section 'Electrical Connection' in the desired configuration (ECS, analog 4–20 mA interface, digital RS485 interface. The combination of RS485 and ECS or RS485 and 4–20 mA interfaces is also possible).

The sensor signal stabilizes itself within few minutes. The sensor is factory pre-calibrated and ready for use. However, if an on-site calibration is necessary, a calibration can be performed via RS485 by means of the ARC View Handheld, personal computer or VISICAL (P/N: 242410).

### NOTE!

*Due to the integrated optoelectronics the sensor has an intrinsic warming which stabilizes itself within a few minutes. If the sensor is not dipped into the liquid, then it displays a reading that is slightly too high.*

## 8. Removal of the sensor

Before removing any sensor from a process, make sure that no process medium can escape during removal.

Unscrew the PG 13.5 thread connection and pull out the sensor. Never turn the sensor at the socket head because by doing so you can loosen the sensor cap and fluid can reach the interior of the sensor.

## 9. Sterilization, autoclaving, CIP procedures

VISIFERM™ DO sensors are carefully designed to withstand the cleaning procedures normally used in biotechnological industries without the need for special precautions. Nevertheless, as is the case with all sensors, it is seen that frequent cleanings or sterilizations with steam or hot and aggressive solutions lead to shortened life span of the sensor. Although for a number of reasons (variations in composition of cleaning agent, temperature gradient and so on) it is not possible to give general life span figures. With steam sterilizations and autoclaving for 30 min at a temperature of 125°C, a life span of more than 50 cycles per exchangeable sensor cap is common.

The contacts must be clean and dry before sensor is connected to the cable.

### ATTENTION!

*The sensor can communicate over digital RS485 interface up to a temperature of 130°C. However, to protect the optoelectronics and the luminophore, no oxygen measurement is performed above 85°C. Therefore, above a temperature of 85°C the sensor cannot be used for control or monitoring.*

## 10. Testing and maintenance

### Calibration

The VISIFERM™ DO sensors provide two kinds of sensor calibration: automatic standard calibration, and product calibration.

The automatic standard calibration is performed over RS485 interface using the ARC View Handheld, VISICAL or the ARC Sensor Configurator freeware (see the sections entitled 'Configuration and monitoring of the sensor' and 'Accessories'). For the product calibration use the ARC View Handheld, or the ARC Sensor Configurator freeware (VISICAL is not supporting the product calibration.)

#### Automatic standard calibration

The VISIFERM™ DO sensor is calibrated at two points: in air and in an oxygen-free environment. During calibration, the sensor examines the correctness and stability of the oxygen and temperature signals. However, the recalibration of the VISIFERM™ DO is very rarely necessary for most applications.

Calibration at point 1:

1. Set an appropriate operator level (Administrator or Specialist).
2. Immerse the sensor into an oxygen-free environment.
3. Let the system equilibrate. Guarantee stable conditions for at least three minutes. For greater measurement accuracy insure that temperature difference between calibration medium and process medium is minimal.
4. Execute the calibration at point 1.
5. If all conditions are met, the sensor confirms the calibration immediately. The calibration curve of the sensor is now defined by the recently calibrated point 1 and the prior calibration at point 2.

Calibration at point 2:

1. Set an appropriate operator level (Administrator or Specialist).
2. Leave the sensor for at least three minutes under stable conditions in ambient air or in oxygen saturated medium.
3. Execute the calibration at point 2.
4. If all conditions are met, the sensor confirms the calibration immediately. The calibration curve of the sensor is now defined by the recent calibration at both points: 1 and 2.

If you experience measurement value shifts as small as a few percent within a period of a few days, suspect damage to the luminophore caused by the measurement or cleaning medium. Should this occur, consider whether more frequent exchange of the sensor cap is acceptable, or whether a different approach to measurement (for example, using a classical OXYFERM sensor) should be taken.

**HINT:** The sensor requires a preheating period of 10 to 15 minutes after it is switched on. Although measurements are possible during this time, you must wait until preheating is complete to be able to calibrate the sensor optimally.

**HINT:** The concept behind VISIFERM™ DO sensors enables calibration and configuration in the lab before use in the process control. Another calibration for the installation in the process setup is not required.

#### **Product calibration.**

The product calibration is an in-process calibration procedure in order to adjust the measurement to specific process conditions, or in case the sensor cannot be removed for the standard calibration.

Product calibration is an additional calibration procedure to a standard calibration. Product calibration adapts the standard calibration curve to the process conditions in force at the time of product calibration. If product calibration is activated, the Visiferm's calibration curve is calculated from the data of last calibration at point 1 and from the data of the product calibration. In order to restore the original standard calibration curve, the product calibration can be cancelled at any time. A new standard calibration cancels a product calibration as well.

Actions:

1. Set an appropriate operator level (Administrator or Specialist).
2. Perform an initial measurement while taking a sample from the process. The data of the initial measurement are stored in the sensor.
3. Perform a laboratory measurement of the sample at the same temperature as it was measured in the process.
4. Assign the laboratory value to the value of the Initial measurement. This new DO value is accepted and instantaneous active, if the difference between initial measurement and laboratory values is not greater than 20%-sat units.

**NOTE!**

*The product calibration is possible for DO values in the range of 2 to 55 %-vol (20 – 550 mbar).*


## Calibration in ECS mode

The calibration in ECS mode takes place like in classical sensors in the measuring device at which the sensor is connected. For this, consult the operation instructions of the measuring device!

**HINT:** The concept of VISIFERM DO sensors enables the pre-calibration using the ARC View Handheld, VISICAL or the ARC Sensor Configurator freeware and then a direct use at classical measuring devices without a new calibration at the measuring device. This is possible, if the sensor was recently calibrated and in the configuration of the sensor's ECS interface an air current of 60 nA was adjusted. The measuring device has set the air current at 25°C on 60 nA correspondingly.

## Changing the sensor cap

The exchange of sensor cap is performed very easily:

Unscrew the old sensor cap from the shaft.

**HINT:** If the sensor cap is mounted very firmly on the shaft, and if you cannot obtain a good grip on the stainless steel with your fingers, a silicone tube between your fingers and metal may supply a better grip.

Examine the small O-ring that seals the sensor cap to the sensor shaft. Exchange the O-ring, if any traces of wear are seen. A replacement O-ring is included with each replacement sensor cap.

Screw the new sensor cap onto the sensor shaft again. Make sure that the gap between the shaft and cap is closed, and that the sealing effect of the O-ring under it is therefore guaranteed. In order to enable traceability (as required by the Federal Drugs Administration) note the serial number of the new sensor cap.

Examine the measurement values of the sensor in air, and if necessary, in an oxygen-free medium. If the measurement values deviate significantly from operated value, perform a calibration.

## Self-diagnosis functions

VISIFERM™ DO sensors provide a self-diagnosis functionality to detect and identify the most common sensor malfunctions. Both interfaces can be used for warning and error messaging. The analog 4–20 mA interface can be configured according to the NAMUR recommendations to indicate an abnormal event. The RS485 interface provides a variety of the indications based on the fault code.

Use the ARC View Handheld and ARC Wi Sensor Adapter for monitoring of the sensor status and troubleshooting. These devices detect and highlight the sensor condition according to the fault code. The ARC View Handheld displays the text messages that are appropriate for the type of malfunction.

The following types of messages are provided by the self-diagnosis function:

- Warning (alarm):
  - DO readings upper/lower range, unstable.
  - Temperature readings upper/lower, unstable.
  - Calibration recommended.
  - Calibration upper/lower, out of range, or unstable.
  - Replace sensor cap.
  - 4–20 mA output out of range, unstable
  - Supply voltage upper/lower range.
- Error (failure):
  - Sensor cap missing (reading failure).
  - Temperature out of range.
  - Temperature sensor defect (reading failure).
  - Supply voltage out of range.

## 11. Disposal

The design of HAMILTON sensors optimally considers environmental compatibility. In accordance with the EC guideline 2002/96/EG HAMILTON sensors that are worn out or no longer required must be sent to a dedicated collection point for electrical and electronic devices, alternatively, must be sent to HAMILTON for disposal. Sensors must not be sent to an unsorted waste disposal point.


## 12. Parts and Accessories

ARC View Handheld with Dock.

ARC Wi Sensor Adapter for wireless communication.

VISICAL.

USB-RS485 Modbus Converter.

Sensor Cable VP8.

Order number	Description
242410	VISICAL
242411	USB-RS485 Modbus converter
242180	ARC View Handheld Package
355194	Demo Cable (1m, open end, with plug power pack)
242170	ARC Wi Sensor Adapter
238999-2394	Sensor cable VP8, 1m
238999-1953	Sensor cable VP8, 3m
238999-2395	Sensor cable VP8, 5m
238999-2396	Sensor cable VP8, 10m
238999-2403	Sensor cable VP8, 15m
238999-2505	Sensor cable VP8, 20m

### 13. Technical data

- VISIFERM™ DO sensor with integrated electronics.  
Functionality includes measurement and self-diagnosis.
- Steam sterilizable, autoclavable and compatible with Cleaning In Place (CIP).
- Measurement principle: oxygen dependent luminescence.
- Exchangeable sensor cap.
- Shaft has a diameter of 12 mm and a PG 13.5 thread.
- Different lengths of shaft available, from 120 mm to 425 mm.
- Operating temperature range: -10 to 130°C; the sensor provides no DO readings above 85°C
- Process pressure range: 0 to 12 bar / 174 psi (relative).
- Range of measurement: 0.05 to 300% air saturation.
- Measured values can be configured by software according to:
  - DO: %-sat, %-vol, mg/l or ppm, ug/l or ppb.
  - Temperature °C.
- Response time at 25°C, from air to nitrogen;  $t_{(98\%)} < 30$  s
- Detection limit: 0.01 %-vol.
- Drift: lower than 0.2 %-vol oxygen per week in air at 30°C and constant conditions.
- Cross sensitivities and resistances:
  - Sensor not disturbed by: carbon dioxide, hydrogen sulfide, sulfur dioxide, ethylene oxide or gamma-sterilization.
  - Sensor resistant to: ethanol, methanol, hydrogen peroxide
  - Sensor not resistant to: chlorine gas and other organic solvents such as chloroform, toluene, acetone.

- Medium affected materials: see certificate.
- Electrical connection: VP8 socket head.
- 7 to 30 VDC operating voltage; max. 1W, 0.6 W continuous.
- One scalable ECS DO interface, or 4–20 mA current interface for DO or temperature signals.
- Digital RS485 Interface (two-wire):
  - Protocol: Modbus RTU; maximal 31 addresses.
  - Baud rate: 4800, 9600, 19200, 38400, 57600, 115200 bd.

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Copyright © 2010 HAMILTON Bonaduz AG, Bonaduz Schweiz. Alle Rechte vorbehalten. Die Reproduktion irgendeines Teils dieses Dokuments in jeder beliebigen Form ist ohne die ausdrückliche schriftliche Zustimmung der HAMILTON Bonaduz AG untersagt.