


# SoftStart™ Pipette Kits

*The best performance and comfort  
at an even better value!*



*3 Racks  
of Free Tips!*



## Adjustable Volume Pipette Kit

- Any **THREE** SoftGrip™ adjustable volume pipettes
- Receive three racks of matching AdvanTip™ pipette tips
- 6 position pipette carousel
- Free shipping

**\$599!**

Value up to  
\$1,100.00!

## Fixed Volume Pipette Kit

- Any **THREE** SoftGrip™ fixed volume pipettes
- Receive three racks of matching AdvanTip™ pipette tips
- 6 position pipette carousel
- Free shipping

**\$399!**

Value up to  
\$790.00!



*10 Racks  
of Free Tips!*



## Multi-Channel Pipette Kit

- Any SoftGrip™ 8-channel or 12-channel pipette
- Any box (10 racks\*) of AdvanTip™ pipette tips
- One pipette stand
- Free shipping

**\$599!**

Value up to  
\$940.00!

**For more information or to place your order,  
call Hamilton at 1-888-525-2123**



**HAMILTON**

Offer ends December 31, 2011 and is only good in the U.S. when purchased directly from Hamilton. Offer cannot be combined with any other Hamilton discounts or offers. Use code 'SSKITS' when placing order.

\*Rack quantity is 96 or 100 tips/rack.

## How to Order

Use the list below to select your pipettes. The list price is shown for reference only. You will receive the special kit pricing by contacting sales at 1(888) 525-2123. For more details on Hamilton pipette tips, go to [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

### SoftGrip™ Adjustable Volume Pipettes

Part No.	Volume	Increments ( $\mu$ L)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume		List Price
				Accuracy within $\pm\%$	Precision within $+\%$	Accuracy within $\pm\%$	Precision within $+\%$	
55019-30	0.2-2 $\mu$ L	0.002	Aqua	8.00	4.00	1.20	0.60	\$285.00
55019-32	1-10 $\mu$ L	0.02	Purple	2.50	1.20	0.80	0.40	\$285.00
55019-34	2.5-25 $\mu$ L	0.02	Forest Green	4.50	1.50	0.80	0.20	\$260.00
55019-36	10-100 $\mu$ L	0.2	Violet	1.80	0.70	0.60	0.15	\$260.00
55019-38	30-300 $\mu$ L	0.2	Brick Red	1.20	0.40	0.40	0.15	\$260.00
55019-40	100 $\mu$ L-1 mL	2.0	Sky Blue	1.60	0.50	0.30	0.12	\$260.00

### SoftGrip™ Fixed Volume Pipettes

Part No.	Volume	Color	At 100% of Pipette Volume		List Price
			Accuracy within $\pm\%$	Precision within $+\%$	
55019-01	5 $\mu$ L	Steel Blue	1.40	0.75	\$170.00
55019-03	10 $\mu$ L	Purple	0.80	0.40	\$170.00
55019-05	25 $\mu$ L	Forest Green	0.80	0.30	\$170.00
55019-07	50 $\mu$ L	Sandstone	0.60	0.30	\$170.00
55019-09	100 $\mu$ L	Violet	0.50	0.20	\$170.00
55019-11	200 $\mu$ L	Mustard	0.40	0.18	\$170.00
55019-13	250 $\mu$ L	Burnt Orange	0.40	0.18	\$170.00
55019-15	300 $\mu$ L	Brick Red	0.40	0.18	\$170.00
55019-17	500 $\mu$ L	Olive Green	0.40	0.18	\$170.00
55019-19	1 mL	Sky Blue	0.30	0.12	\$170.00

### SoftGrip™ Multi-Channel Pipettes, 8-Channel

Part No.	Volume	Increments ( $\mu$ L)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume		List Price
				Accuracy within $\pm\%$	Precision within $+\%$	Accuracy within $\pm\%$	Precision within $+\%$	
55019-50	5-50 $\mu$ L x 8	0.02	Sandstone	2.00	2.00	0.80	0.35	\$680.00
55019-52	30-300 $\mu$ L x 8	0.2	Brick Red	1.20	0.70	0.50	0.30	\$680.00

### SoftGrip™ Multi-Channel Pipettes, 12-Channel

55019-54	5-50 $\mu$ L x 12	0.02	Sandstone	5.00	2.00	1.00	0.50	\$800.00
55019-56	30-300 $\mu$ L x 12	0.2	Brick Red	2.40	1.20	1.00	0.50	\$800.00

### SoftGrip™ Pipette Carousel

Keeps pipettes safe and within easy reach. Constructed of heavy, solvent-resistant plastic.

Description	Part No.	List Price
6 Position Pipette Carousel	53578-01	\$80.00



### SoftGrip™ Multi-Channel Pipette Stand

Keeps pipette in easy reach. Constructed of a heavy, non-skid, metallic casting.

Description	Part No.	List Price
Holds One Multi-Channel Pipette	34195	\$77.30



**For more information or to place your order,  
call Hamilton at 1-888-525-2123**

Offer ends December 31, 2011 and is only good in the U.S. when purchased directly from Hamilton. Offer cannot be combined with any other Hamilton discounts or offers. Use code 'SSKITS' when placing order.



Hamilton Company  
4970 Energy Way  
Reno, Nevada 89502 USA  
Toll-Free: 800-648-5950  
Telephone: +1-775-858-3000  
Fax: +1-775-856-7259  
Email: [sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

Hamilton Bonaduz AG  
Via Crusch 8  
CH-7402  
Bonaduz/Switzerland  
Telephone: +41-(0)81-660-6060  
Fax: +41-(0)81-660-6070  
Email: [contact@hamilton.ch](mailto:contact@hamilton.ch)

[www.hamiltoncompany.com](http://www.hamiltoncompany.com)