

WOW!

SofTouch™ Electronic Pipettes
starting at

\$299



Simplify tasks with seven easy-to-program pipetting modes and five aspirate and dispense speed settings while minimizing hand fatigue with soft action and lightweight pipette design.

- Order your choice of any SofTouch electronic single channel or multi-channel pipette
- Choose one **FREE** box (10 racks) of AdvanTip™ Pipette Tips
- **FREE** shipping

For more information or to place your order, call Hamilton at 1-888-525-2123



HAMILTON

Offer ends December 31, 2009 and is only good in the U.S. when purchased directly from Hamilton. Offer cannot be combined with any other Hamilton discounts or offers. Use code 'WOW' when placing order.

Selecting Hamilton SofTouch Pipettes

Use the list below to select your pipettes. The list price is shown for reference only. You will be receiving the special pricing. For more details on SofTouch pipettes, go to www.hamiltoncompany.com

SofTouch Single Channel Electronic Pipettes, 0.2 µL - 5 mL

Part No.	Volume	Increments (µL)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume		List Price	Your Price
				Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%		
8707-01US	0.2-10 µL	0.1	Purple	2.50	1.50	0.90	0.50	\$530.00	\$299.00
8707-02US	5-100 µL	1.0	Violet	2.00	1.00	0.40	0.15	\$530.00	\$299.00
8707-03US	50-1000 µL	5.0	Sky Blue	1.50	1.50	0.40	0.15	\$530.00	\$299.00
8707-04US	100-5000 µL	50.0	Graphite	0.80	0.30	0.50	0.15	\$530.00	\$299.00

SofTouch Multi-Channel Electronic Pipettes, 8-Channel, 0.2 µL - 1200 µL

Part No.	Volume	Increments (µL)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume		List Price	Your Price
				Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%		
8707-05US	0.2-10 µL	0.1	Purple	4.00	4.00	0.90	0.50	\$998.00	\$549.00
8707-06US	5-100 µL	1.0	Violet	2.50	1.50	0.50	0.20	\$947.00	\$549.00
8707-07US	25-250 µL	5.0	Burnt Orange	1.50	1.00	0.40	0.10	\$947.00	\$549.00
8707-08US	50-1200 µL	10.0	Teal	4.00	0.80	0.50	0.15	\$998.00	\$549.00

SofTouch Multi-Channel Electronic Pipettes, 12-Channel, 0.2 µL - 1200 µL

Part No.	Volume	Increments (µL)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume		List Price	Your Price
				Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%		
8707-09US	0.2-10 µL	0.1	Purple	4.00	4.00	0.90	0.50	\$1,109.00	\$649.00
8707-11US	25-250 µL	5.0	Burnt Orange	1.50	1.00	0.40	0.10	\$1,031.00	\$649.00
8707-12US	50-1200 µL	10.0	Teal	4.00	0.80	0.80	0.15	\$1,109.00	\$649.00

Selecting Hamilton AdvanTip Pipette Tips

For details on Hamilton pipette tips, go to www.hamiltoncompany.com

Identify the pipette tips that will best suit your needs. You will receive one free box of 10 racks with each pipette. Hamilton pipette tips are universal-fit.

Hamilton offers two high quality tip families: AdvanTip Precision tips for routine, high-throughput work and AdvanTip LT (low retention) tips for maximizing sample recovery. Multiple formats and packaging options are available. All AdvanTip products are certified to be free of RNase, DNase, and pyrogens. Hamilton AdvanTip Pipette Tips provide outstanding accuracy and precision, guaranteed.



**Buy direct and save over 50%,
call Hamilton at 1-888-525-2123**

Offer ends December 31, 2009 and is only good in the U.S. when purchased directly from Hamilton. Offer cannot be combined with any other Hamilton discounts or offers. Use code 'WOW' when placing order.



Hamilton Company
4970 Energy Way
Reno, Nevada 89502 USA
Toll-Free: 800-648-5950
Telephone: +1-775-858-3000
Fax: +1-775-856-7259
Email: sales@hamiltoncompany.com

Hamilton Bonaduz AG
Via Crusch 8
CH-7402
Bonaduz/Switzerland
Telephone: 00800-660-660-60
Fax: +41-(0)81-660-60-70
Email: marketing@hamilton.ch

www.hamiltoncompany.com