

Microlab 600 Universal Valve - Left

Specification Sheet (Part/REF # 60676-01)



The Microlab Universal Valves are used on the Single Syringe Dispenser, Dual Syringe Diluter, Dual Syringe Dispenser and Continuous Dispenser. Adjust the Cross Tube and Valve Plugs to achieve the desired application.

Product Specifications

Valve Model Description	Microlab 600 Universal Valve Left
Fluid Path	PTFE/CTFE
Nominal Port Diameter	0.059 inch (1.5 mm)
Housing Material	PEEK
Actuated	All Microlab 600 instruments
Break Away Torque	65 in. oz. (455 mNm)
Panel Mounting	No
Threaded Connection	1/4-28
Operating Temp High	38 °C (100 °F)
Operating Temp Low	-12 °C (10 °F)
Operating Pressure	0 to 100 psig (0 to 6.9 bar)
Storage Temp High	38 °C (100 °F)
Storage Temp Low	-18 °C (0 °F)

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/valves/60676-01>

Spec. Version B