

# Microlab 600 Universal Valve - Left

Specification Sheet (Part/REF # 60676-01)



The Microlab Universal Valves are used on the Single Syringe Dispenser, Dual Syringe Diluter, Dual Syringe Dispenser and Continuous Dispenser. Adjust the Cross Tube and Valve Plugs to achieve the desired application.

## Product Specifications

<b>Valve Model Description</b>	Microlab 600 Universal Valve Left
<b>Fluid Path</b>	PTFE/CTFE
<b>Nominal Port Diameter</b>	0.059 inch (1.5 mm)
<b>Housing Material</b>	PEEK
<b>Actuated</b>	All Microlab 600 instruments
<b>Break Away Torque</b>	65 in oz. (455 mNm)
<b>Panel Mounting</b>	No
<b>Threaded Connection</b>	¼-28
<b>Operating Temp High</b>	38 °C (100 °F)
<b>Operating Temp Low</b>	-12 °C (10 °F)
<b>Operating Pressure</b>	0 to 100 psig (0 to 6.9 bar)
<b>Storage Temp High</b>	38 °C (100 °F)
<b>Storage Temp Low</b>	-18 °C (0 °F)

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/valves/60676-01>

**Spec. Version B**