

Hamilton MIC-Tubes and V-Caps - Only available in Europe

Specification Sheet (Part/REF # 10110770)



Hamilton MIC Tubes and V-Caps are used with the MIC4Hamilton (Magnetic Induction Cycler). During use the MIC4Hamilton is fully integrated on a Hamilton automated liquid handler. Pre-assembled V-Caps provides workflow efficiency with reduced tube reload time.

Product Specifications

Compatibility	MIC4Hamilton (Magnetic Induction Cycler)
Volume Format	0.1 mL
Volume supported	5 - 30 µL
Sterile	Non-Sterile
Tube Material	Polypropylene
V-CAP Material	Polypropylene
Rack & Lid Material	Polybutylene Terephthalate (PBT, upcycled material)
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, PCR inhibitors and ATP
IVD-R Compliant	Yes
Tubes/Rack	48 pieces
V-CAPS/Rack	48 pieces
Racks/Case	20 pieces
Case size	960 Tubes and V-CAPS

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/10110770>

Spec. Version B