

1.2 mL Deep Well Plate With Barcode Label - Only Available in the U.S.

Specification Sheet (Part/REF # 235655)



1.2 mL Deep Well Plate; Polystyrene with Barcode Label

Product Specifications

Maximum Well Capacity	1.2 mL
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, ATP, and PCR inhibitors
Material	Polystyrene - PS
Recycling Code	6
SLAS ANSI	Yes
Packaging	4 plates per bag
Barcode	Yes
Compatible Carrier	SLAS/ANSI carrier
Quantity	32 plates per case

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/235655>

Spec. Version B