

Sealing Mat for 2.2 mL Deep Well Plate

Specification Sheet (Part/REF # 6474-01)



Sealing Mat for 2.2 mL Deep Well Plates

Product Specifications

Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, ATP, and PCR inhibitors
Material	Ethylene vinylacetate (EVA)
Packaging	10 mats per bag
Barcode	No
Quantity Per Case	50 mats per case
Compatible Carrier	N/A

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/6474-01>

Spec. Version B