

AutoLys Tubes (500 tubes)

Specification Sheet (Part/REF # 235478)



Product Specifications

Purity Level	Sterile Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, PCR inhibitors, and ATP
Material	Polypropylene - PP
Recycling Code	5
SLAS ANSI	N/A
Barcode	Yes
Packaging	100 tubes per bag
Quantity Per Case	500 tubes per case

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/235478>

Spec. Version B