

# 1.5 mL FlipTube Biological purity

Specification Sheet (Part/REF # 235692)



## Product Specifications

<b>Purity Level</b>	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
<b>Material</b>	Polypropylene - PP
<b>Recycling Code</b>	5
<b>Packaging</b>	500 tubes per bag
<b>Quantity Per Case</b>	2000 tubes per bag

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/235692>

**Spec. Version B**