

1.5 mL FlipTube Biological Purityplus Sterile

Specification Sheet (Part/REF # 235693)



Product Specifications

Purity Level	Sterile Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, PCR inhibitors, and ATP
Material	Polypropylene - PP
Recycling Code	5
Packaging	100 tubes per bag
Quantity Per Case	2000 tubes per case

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/235693>

Spec. Version B