

60 mL Reagent Reservoir, Self- Standing with Lid

Specification Sheet (Part/REF # 56694-01)



60 mL Reagent Reservoir, Self-Standing with Lid, No Barcode

Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produced in a Class 8 Clean Room
Purity Level	Free of DNase, RNase, Human DNA, and PCR inhibitors
Recycling Code	5
Quantity Per Case	Case of 28 reservoirs
Packaging	Individually wrapped
SLSA ANSI	No
Barcode	No
Compatible Carrier	53646-01 60 mL Reagent Reservoir Carrier - 5 reservoir positions

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/56694-01>

Spec. Version B