

# 60 mL Black Reagent Reservoir Self- Standing with Lid

Specification Sheet (Part/REF # 56694-02)



60 mL Black Reagent Reservoir Self-Standing with Lid, No Barcode

## Product Specifications

|                       |  |
|-----------------------|--|
| Material              | Polypropylene - PP   |
| Production Compliance | Product produced in a Class 8 Clean Room                         |
| Purity Level          | Free of DNase, RNase, Human DNA, and PCR inhibitors              |
| Recycling Code        | 5  |
| Quantity              | Case of 28 reservoirs  |
| Packaging             | Individually wrapped   |
| SLAS ANSI             | No   |
| Barcode               | No   |
| Compatible Carrier    | 53646-01 60 mL Reagent Reservoir Carrier - 5 reservoir positions |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/56694-02>

Spec. Version B