

# 60 mL Reagent Reservoir Self- Standing

Specification Sheet (Part/REF # 194051)



60 mL Reagent Reservoir, Self-Standing without Lid, Barcode Label included

## Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produced in a Class 8 Clean Room
Purity Level	Free of DNase, RNase, Human DNA, and PCR inhibitors
Recycling Code	5
Quantity Per Case	Case of 28 reservoirs
Packaging	Individually wrapped
SLAS ANSI	No
Barcode	Yes
Compatible Carrier	53646-01 60 mL Reagent Reservoir Carrier - 5 reservoir positions

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/194051>

Spec. Version B