

200 mL Reagent Reservoir, Self- Standing with Lid

Specification Sheet (Part/REF # 56695-01)



200 mL Reagent Reservoir, Self-Standing with Lid, No Barcode Label

Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produced in a Class 8 Clean Room
Purity Level	Free of DNase, RNase, Human DNA, and PCR inhibitors
Recycling Code	5
Quantity Per Case	Case of 10 reagent reservoirs
Packaging	Individually wrapped
SLAS ANSI	No
Barcode	No
Compatible Carrier	53645-01 200 mL Reagent Reservoir Carrier - 4 reservoir positions
Compatible Aluminum Lid	p/n 6603425-01

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/56695-01>

Spec. Version B