

# 200 mL Black Reagent Reservoir, Self-Standing with Lid

Specification Sheet (Part/REF # 56695-02)



200 mL Black Reagent Reservoir, Self-Standing with Lid, No Barcode Label

## Product Specifications

Material	Polypropylene - PP
Production Compliance	Product produce in a Class 8 Clean Room
Purity Level	Free of DNase, RNase, Human DNA, and PCR inhibitors
Recycling code	5
Quantity Per Case	Case of 10 reagent reservoirs
Packaging	Individually wrapped
Barcode	No
SLAS ANSI	No
Compatible Carrier	53645-01 200 mL Reagent Reservoir Carrier - 4 reservoir positions

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/56695-02>

Spec. Version B