


2.2 mL Deep Well Plate with Barcode Label

Specification Sheet (Part/REF # 235656)



2.2 mL Deep Well Plate with Barcode Label

Product Specifications

Maximum Well Capacity	2.2 mL
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, ATP, and PCR inhibitors
Material	Polypropylene - PP
Recycling Code	5
Autoclavable	Yes
Well Bottom	Round
Well Format	96
SLAS ANSI	Yes
Barcode	Yes
Compatible Carrier	SLAS/ANSI Carrier
Quantity	32 plates per case

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/automated-liquid-handling/other-robotics/part-ref-235656>