

# Microlab 600 Single Syringe Pump, PC Controlled

Specification Sheet (Part/REF # ML630)



Program methods and deploy commands from your PC via Ethernet or RS-232 communication for unparalleled control of your process.

## Product Specifications

<b>Power Rating</b>	24 VDC, 2.5A
<b>Power Requirements</b>	100-240V 1.5A max 50-60 Hz
<b>Program Memory</b>	One program stored in memory while power is on
<b>Syringe Drive Mechanism</b>	Stepper motor driven high precision lead screw with encoder
<b>Controller</b>	PC Controlled
<b>Number of Syringe Drives</b>	One
<b>Probe</b>	No Probe

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml630>