

Microlab 600 Dual Syringe Pump, PC Controlled

Specification Sheet (Part/REF # ML635)



Includes the dual syringe pump, two valves, cross tube assembly, two valve plugs, power supply, country-specific power cord, Custom Programmer Kit, and your choice of two syringes. If syringe volumes are not selected at the time of order, default 1 mL syringes will be included automatically.

Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Program Memory	One program stored in memory while power is on
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder
Controller	PC Controlled
Number of Syringe Drives	Two
Probe	No Probe
Fluid Path	Borosilicate, PTFE, CTFE, PFA
Dimensions	7 x 5.5 x 10.5 inch (177.8 x 139.7 x 266.7 mm)
Weight	13 lbs

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml635>

Spec. Version B