

Microlab 600 Advanced Dual Syringe Dispenser

Specification Sheet (Part/REF # ML625-DIS)



The Dual Syringe Dispenser has two syringe drives for applications such as Animal Dosing that require the use of more than one liquid at the same time. The two liquids can be dosed from the same instrument saving setup and dosing time.

Product Specifications

| | |
|--------------------------|---|
| Power Rating | 24 VDC, 2.5A |
| Power Requirements | 100-240V 1.5A max 50-60 Hz |
| Program Memory | 2 GB |
| Valve Setup | Dual Syringe Dispenser |
| Syringe Drive Mechanism | Stepper motor driven high precision lead screw with encoder |
| Controller | Advanced |
| Number of Syringe Drives | Two |
| Probe | Dual Push Button Hand Probe |
| Syringes Included | Includes two user selected syringes |
| Fluid Path Materials | FEP, CTFE, PTFE, Borosilicate Glass |

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml625-dis>

Spec. Version B