

Microlab 600 Advanced Continuous Dispenser

Specification Sheet (Part/REF # ML625-CNT)



The Continuous Dispenser is designed for time critical applications requiring repetitive dispensing. One syringe fills while the other syringe dispenses so the instrument is always ready. Perfect for integration into semi-automated or fully automated manufacturing lines.

Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Program Memory	2 GB
Valve Setup	Continuous Dispenser
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder
Controller	Advanced
Number of Syringe Drives	Two
Probe	Concorde CT Hand Probe
Syringes Included	Includes two user selected syringes
Fluid Path Materials	FEP, CTFE, PTFE, Borosilicate Glass

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml625-cnt>

Spec. Version B