

Analytical Guard Column, PRP-1

Specification Sheet (Part/REF # 79286)



Protect your PRP-1 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
USP Classification	USP L21
Packing Material Type	PSDVB
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79286>