

PRP-X100 Anion Exchange Analytical Guard Column

Specification Sheet (Part/REF # 79287)



Protect your PRP-X100 column.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X100
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
Packing Material Type	PSDVB/Trimethylammonium
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79287>