


# Analytical Guard Column, PRP-X300

Specification Sheet (Part/REF # 79289)



Organic acids and alcohols.

## Product Specifications

<b>Separation Mode</b>	Ion Exclusion
<b>Stationary Phase</b>	PRP-X300
<b>Format</b>	Guard Column
<b>Hardware Inner Diameter</b>	2.1 mm
<b>Hardware Length</b>	33 mm
<b>Particle Size</b>	12-20 µm
<b>Packing Material Type</b>	PSDVB/Sulfonic Acid
<b>Guard Column Size</b>	Analytical
<b>Hardware Material</b>	Stainless Steel
<b>Max Pressure</b>	5,000 psi
<b>Compound Classification</b>	Organic Acids/Alcohols
<b>UOM</b>	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79289>