

Analytical Guard Column, HxSil C8

Specification Sheet (Part/REF # 79298)



Polycyclic aromatic hydrocarbons, small organic molecules. Shorter retention than C18.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	HxSil C8
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
Packing Material Type	Silica
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Peptides/Proteins
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79298>