

Analytical Guard Column, PRP-3

Specification Sheet (Part/REF # 79295)



Organic compounds: large molecules (> 2,000 mw), peptides, proteins, protein digests, protected and deprotected oligonucleotides, nucleic acids.

Product Specifications

| | |
|--------------------------------|-----------------|
| Separation Mode | Reversed-Phase |
| Stationary Phase | PRP-3 |
| Format | Guard Column |
| Hardware Inner Diameter | 2.1 mm |
| Hardware Length | 33 mm |
| Particle Size | 10 µm |
| Packing Material Type | PSDVB |
| Guard Column Size | Analytical |
| Hardware Material | Stainless Steel |
| Max Pressure | 5,000 psi |
| Compound Classification | Organics |
| UOM | 1 EACH |

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79295>