

# Analytical Guard Column, PRP-h5

Specification Sheet (Part/REF # 79296)



Organic compounds: large molecules (> 2,000 mw), halogenated compounds, pharmaceuticals, protein digests, tryptic digests, proteomics.

## Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-h5
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
Packing Material Type	Fluorinated PSDVB
Particle Surface Area (m <sup>2</sup> /cm <sup>3</sup> )	277
Pore Volume (cm <sup>3</sup> /g)	0.526
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Peptides/Proteins
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79296>