

Analytical Guard Column, PRP-X400

Specification Sheet (Part/REF # 79290)



Glyphosate and its metabolites in drinking water, inorganic and organic cations.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X400
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	7 µm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Glyphosate
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79290>