

Semiprep/Preparative Guard Column, PRP-1

Specification Sheet (Part/REF # 79912)



Organic compounds: small molecules (< 2,000 mw), oligonucleotides, pharmaceuticals, steroids, nucleic acids, vitamins, herbicides.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Column
Hardware Inner Diameter	4.6 mm
Hardware Length	33 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organics
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79912>