

# Semiprep/Preparative Guard Column, PRP-C18

Specification Sheet (Part/REF # 79918)



Organic compounds: Oligonucleotides, small molecules (< 2,000 mw), pharmaceuticals, steroids, halides, vitamins, amino acids, herbicides. High pH applications.

## Product Specifications

<b>Separation Mode</b>	Reversed-Phase
<b>Stationary Phase</b>	PRP-C18
<b>Format</b>	Guard Column
<b>Hardware Inner Diameter</b>	4.6 mm
<b>Hardware Length</b>	33 mm
<b>Particle Size</b>	12-20 $\mu$ m
<b>Packing Material Type</b>	PS-DVB functionalized with C18
<b>Guard Column Size</b>	Semiprep/Preparative
<b>Hardware Material</b>	Stainless Steel
<b>Max Pressure</b>	5,000 psi
<b>Compound Classification</b>	Organics
<b>UOM</b>	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79918>