

Semiprep/Preparative Guard Column, PRP-h5

Specification Sheet (Part/REF # 79921)



Organic compounds: large molecules (> 2,000 mw), halogenated compounds, pharmaceuticals, protein digests, tryptic digests, proteomics.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-h5
Format	Guard Column
Hardware Inner Diameter	4.6 mm
Hardware Length	33 mm
Particle Size	12-20 µm
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Peptides/Proteins
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79921>