

# PRP-X100 Anion Exchange Semiprep/Preparative Guard Column

Specification Sheet (Part/REF # 79913)



High capacity ion exchange chromatography of organic and inorganic anions, organic acids, organic and inorganic arsenic species, halides, nucleotides.

## Product Specifications

<b>Separation Mode</b>	Anion Exchange
<b>Stationary Phase</b>	PRP-X100
<b>Format</b>	Guard Column
<b>Hardware Inner Diameter</b>	4.6 mm
<b>Hardware Length</b>	33 mm
<b>Particle Size</b>	12-20 $\mu\text{m}$
<b>Packing Material Type</b>	PSDVB/Trimethylammonium
<b>Guard Column Size</b>	Semiprep/Preparative
<b>Hardware Material</b>	Stainless Steel
<b>Max Pressure</b>	5,000 psi
<b>Compound Classification</b>	Organic/Inorganic Anions
<b>UOM</b>	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79913>