

Semiprep/Preparative Guard Column, PRP-X300

Specification Sheet (Part/REF # 79915)



Organic acids and alcohols.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Guard Column
Hardware Inner Diameter	4.6 mm
Hardware Length	33 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organic Acids/Alcohols
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79915>