


# Septa, Low Bleed High Temp., Plug style for Shimadzu instruments (50/PK)

Specification Sheet (Part/REF # 75837)



Low Bleed High Temperature Septa are designed for the lowest bleed possible at temperatures up to 400 °C. They can be used in combination with Mass Spectrometers or in applications where septa bleed interferes with trace analysis.

## Product Specifications

Diameter	Plug style for Shimadzu instruments
Maximum Temperature	300°C / 572°F
Thickness	7.00 mm (0.28 inch)
Sales Quantity	50/pack

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/gc-septa/75837>

Spec. Version B