

Septa, 3-Layer, 19.0 mm DIA (12/PK)

Specification Sheet (Part/REF # 76020)



Three Layer Septa is made of preconditioned silicone rubber and is 3.81 mm (0.150 in) thick. The average life of these septa is up to 100 injections with a 26s gauge needle. The maximum operating temperature is 200° C.

Product Specifications

| | |
|---------------------|---------------------------|
| Diameter | 19 mm (0.748 inch) |
| Maximum Temperature | 200°C / 392°F |
| Material | Laminated Silicone Rubber |
| Thickness | 3.81 mm (0.150 inch) |
| Sales Quantity | 12 septa per pack |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/gc-septa/76020>

Spec. Version B