

Carbonate Form Semiprep/Preparative Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79866)



Protect your HC-75 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X500
Format	Guard Cartridge Kit
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	400 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Carbohydrates
UOM	1 holder, 1 cartridge

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79866>