

# Analytical Guard Cartridge Holder, Stainless Steel

Specification Sheet (Part/REF # 32908)



High exchange capacity anion exchange column for carbohydrate separation: Carbohydrates, polysaccharides, sugar oligomers up to DP8.

## Product Specifications

Format	Guard Cartridge Holder
Guard Column Size	Analytical
Hardware Material	Stainless Steel
UOM	1 EACH

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/32908>

Spec. Version B