

PRP-X100 Anion Exchange Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79448)



Protect your PRP-X100 column.

Product Specifications

| | |
|--------------------------------|---|
| Separation Mode | Anion Exchange |
| Stationary Phase | PRP-X100 |
| Format | Guard Cartridge Kit |
| Hardware Inner Diameter | 2.1 mm |
| Hardware Length | 20 mm |
| Particle Size | 10 µm |
| Packing Material Type | PSDVB/Trimethylammonium |
| Guard Column Size | Analytical |
| Hardware Material | Stainless Steel |
| Max Pressure | 5,000 psi |
| Mobile Phase Limits | pH 1-13. 0-100% aqueous, organic modifier |
| Compound Classification | Organic/Inorganic Anions |
| UOM | 1 holder, 2 cartridges |

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79448>