

Analytical Guard Cartridge Holder, PEEK

Specification Sheet (Part/REF # 79477)



High exchange capacity anion exchange column for carbohydrate separation: Carbohydrates, polysaccharides, sugar oligomers up to DP8.

Product Specifications

Format	Guard Cartridge Holder
Guard Column Size	Analytical
Hardware Material	PEEK
UOM	1 EACH

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79477>