

PRP-X100 Anion Exchange HPLC Column, 2.1 x 250 mm, 5 μm

Specification Sheet (Part/REF # 79190)



High capacity ion exchange chromatography of organic and inorganic anions, organic acids, organic and inorganic arsenic species, halides, nucleotides.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X100
Format	HPLC Column
Hardware Inner Diameter	2.1 mm
Hardware Length	250 mm
Particle Size	5 μm
Packing Material Type	PSDVB/Trimethylammonium
Pore Size	100 Å
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Exchange Capacity	0.19 meq/gm
Exchange Capacity/Column	0.076 meq/column
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Buffer Strength	0.0-0.5N
Temp Limits	pH 1-7.9 (5-60°C); pH 8-13 (5-30°C)
Restoration	Pump approximately 50 mL of methanol with 1% 6 N nitric acid.

Compound Classification Organic/Inorganic Anions

UOM 1 EACH

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79190>