

Eluent Concentrate, PRP-X100 Anion Exchange [makes 1.0 L of mobile phase, (60 mL bottle, 1 each)]

Specification Sheet (Part/REF # 79325)



Dilute to 1.0 L with deionized water

Product Specifications

Separation Mode	Anion Exchange
Format	Eluent Concentrate
Sales Quantity	1 x 60 mL bottle (dilute with 1.0 L of deionized water)
Composition	4.0 mM p-Hydroxybenzoic Acid / 2.5% Methanol

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79325>

Spec. Version B