

Eluent Concentrate, PRP-X100 Anion Exchange [each bottle makes 1.0 L of mobile phase] (60 mL bottle, 6 each)

Specification Sheet (Part/REF # 79335)



Quickly prepare mobile phase for the PRP-X100.

Product Specifications

Separation Mode	Anion Exchange
Format	Eluent Concentrate
Sales Quantity	6 x 60 mL bottles (dilute each bottle with 1.0 L of deionized water)
Composition	4.0 mM p-Hydroxybenzoic Acid / 2.5% Methanol

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79335>

Spec. Version B