

PRP-X110 Anion Exchange Analytical Guard Column Replacement Cartridges (5/pk), Stainless Steel

Specification Sheet (Part/REF # 79728)



High sensitivity detection ion exchange chromatography of organic and inorganic anions, organic acids, organic and inorganic arsenic species, halides, nucleotides.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X110
Format	Replacement Guard Cartridge
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Trimethylammonium
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic/Inorganic Anions
UOM	5/pk

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79728>