

PRP-X600 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79361)



Protect your PRP-X600 column.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X600
Format	Guard Cartridge Kit
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Poly (dimethyl amido-propyl methacrylamide)
Particle Surface Area (m ² /cm ³)	1.4
Pore Volume (cm ³ /g)	0.0046
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Labeled/Unlabeled DNA
Frit Material	Titanium
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79361>

Spec. Version B