

Analytical HPLC Column [Custom Order] (1.0, 2.1, 4.1, 4.6 mm ID)

Specification Sheet (Part/REF # 79641)



Organic and inorganic anions, and arsenic speciation.

Product Specifications

Separation Mode	Reversed-Phase Anion Exchange Cation Exchange Ion Exclusion
Stationary Phase	(Choose a phase)
Format	HPLC Column
Hardware Inner Diameter	1.0 mm 2.1 mm 4.1 mm 4.6 mm
Hardware Length	(Choose a length)
Particle Size	(Choose a particle size)
USP Classification	USP L1 USP L7 USP L47 USP L17 USP L19 USP L21 USP L22 USP L23 USP L34
Hardware Material	(Choose hardware material)

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79641>

Spec. Version B