

# Preparative HPLC Column #1, [Custom Order] (21.2, 30 mm ID)

Specification Sheet (Part/REF # 79643)



Inherent reversed-phase characteristics make the PRP-1 a rugged, multi-purpose column.

## Product Specifications

Separation Mode	Reversed-Phase Anion Exchange Cation Exchange Ion Exclusion
Stationary Phase	(Choose a phase)
Format	HPLC Column
Hardware Inner Diameter	21.2 mm 30 mm
Hardware Length	(Choose a length)
Particle Size	(Choose a particle size)
Hardware Material	Stainless Steel

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79643>

Spec. Version B