

RCX-30 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79371)



Increased exchange capacity anion exchange columns for better separation of carbohydrate mixtures: Complex mono- and disaccharide separation.

Product Specifications

| | |
|--------------------------------|---|
| Separation Mode | Anion Exchange |
| Stationary Phase | RCX-30 |
| Format | Guard Cartridge Kit |
| Hardware Inner Diameter | 3.0 mm |
| Hardware Length | 8 mm |
| Particle Size | 12-20 µm |
| Packing Material Type | PSDVB/Trimethylammonium |
| Guard Column Size | Analytical |
| Hardware Material | PEEK |
| Max Pressure | 5,000 psi |
| Mobile Phase Limits | pH 1-13. 0-100% aqueous, organic modifier |
| Compound Classification | Complex Carbohydrates |
| Frit Material | Titanium |
| UOM | 1 holder, 2 cartridges |

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79371>