

RCX-30 Analytical Guard Column Replacement Cartridges (5/pk), PEEK

Specification Sheet (Part/REF # 79372)



Increased exchange capacity anion exchange columns for better separation of carbohydrate mixtures: Complex mono- and disaccharide separation.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	RCX-30
Format	Replacement Guard Cartridge
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 μ m
Packing Material Type	PSDVB/Trimethylammonium
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Complex Carbohydrates
Frit Material	Titanium
UOM	5/pk

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79372>